

Lab Desiccator

Replacement Cartridge

Advantages

- Clear, perforated see-thru Mylar top
- Long-wearing, permeable felt bag
- Cartridges to fit most popular lab desiccators
- Cartridges available with multiple types of desiccant charge: Silica Gel (indicating and white) and Molecular Sieve

Applications

- Lab desiccator jars

Lab Desiccator Replacement Cartridges

The Lab Desiccator Replacement Cartridge was designed to eliminate unwanted mess while changing the desiccant in your lab desiccator jars. It is a replacement cartridge with NatraSorb® indicating desiccant beads that tell you when they are working.



Easy, Convenient, & Ready for Use

Traditional desiccants used in lab desiccator jars are made of loose powder creating a mess and requiring a significant amount of clean up. Unlike conventional desiccants, Lab Desiccator Cartridges can be easily removed and replaced requiring no clean up because desiccant beads are enclosed in a heat sealed felt bag. Indicating silica gel used in lab desiccator cartridges, changes from blue to pink to let you know when your cartridge needs to be replaced.



IMPORTANT NOTICE: Manufacturer/Seller (collectively "Seller") disclaims any warranty for the product sample ("Sample") which Sample is provided "AS IS" without warranty of any kind, either express or implied, including without limitation, the implied warranties of merchantability, fitness for a particular purpose or non-infringement. Seller does not guarantee the accuracy or completeness of any statements made regarding the Sample. User assumes all risk and liability in using the Sample. SELLER SHALL NOT BE LIABLE IN TORT OR CONTRACT FOR ANY LOSS OR DAMAGE, DIRECT, INCIDENTAL, OR CONSEQUENTIAL, ARISING OUT OF THE USE OF OR INABILITY TO USE THE SAMPLE. USER agrees to defend and indemnify Seller from all claims, liabilities, and costs (including attorneys' fees) related to User's use of the Sample. User is prohibited from re-selling Sample.

© 2009 Multisorb Technologies



MULTISORB
Technologies
www.multisorb.com

In the Americas and Asia contact:
Multisorb Technologies
325 Harlem Road
Buffalo, New York 14224 USA
+1 716.824.8900
info@multisorb.com

In Europe and Africa contact:
Multisorb Technologies
Stafford Park 6
Telford
Shropshire
TF3 3AT
United Kingdom
+44 (0) 1952.236.300